

Polar–Non Polar

Purpose

To demonstrate the change in solubility of alcohols in water as the chain length increases.

Materials

ethanol	4 Petri dishes
n-butanol	wash bottle (water)
t-butanol	4 disposable pipets
heptanol	octanol/water bottle (in prep room, prepared
non-polar dye (Sudan)	and colored)

Procedure

1. Place 4 Petri dishes on the overhead.
2. Cover the bottom of all the Petri dishes with water.
3. Add a few drops of each alcohol separately to a Petri dish.
4. Observe the solubility.
5. Ask students to predict octanol and water solubility. Show them the prepared bottle of octanol/water from the prep room.

Additional Information

1. As alcohols become more nonpolar, solubility decreases in water.
2. Sudan dye works well on skin – handle it carefully. If you “dye” yourself a little acetone will reduce the intensity.

Disposal

Solutions containing octanol, n-butanol, t-butanol and heptanol should be placed in a properly labeled waste container.